




ViewⁿVivo
• **See Life Connected**

Handheld, miniaturized, non-destructive, slide-free, confocal endomicroscopy for live preclinical and translational imaging at the single cell level

Optiscan^o



Features

ViewnVivo® by Optiscan is the most advanced miniaturised confocal endomicroscope in the world.

ViewnVivo® allows you to push the boundaries of your research with maximum flexibility, viewing live tissue at any angle with sub-micron resolution.



Miniaturised Hand-held Flexible Probes

Ability to image difficult to access tissue from any angle.




The unique flexible, hand-held ViewnVivo® probe allows stable cellular imaging upon contact of tissue.



Optical Sectioning

Thin optical sections enable 3D reconstruction of tissue structure in remarkable detail.




Optical sectioning of lymphatic vessel captured with ViewnVivo®. Tissue was labeled in vivo with endothelial hyaluronan receptor I (Alexa Fluor 488 conjugate).



Live Pre-clinical and Translational Imaging

4mm diameter probes are designed for non invasive procedures.



Mouse intestine stained with Acriflavine

Mouse liver blood vessels stained with Tomato Lectin.



Probe Choices

Various form factors to suit different research needs.




Different size form factors to suit different research needs.

Ability to continuously image from any angle.




Virtual Histology

ViewnVivo® provides a cutting-edge tool for en-face real-time histological data acquisition.



Endomicroscopy




Histology

Breast Cancer


Fibrous Tissue

Normal Breast


Applications




Transgenic mouse expressing YFP in the cytoplasm of cardiac myocytes of the heart atrium.




Mouse kidney stained with Acriflavine.




Mouse blood vessels stained with i.v. FITC dextran (plasma stain) and acriflavine to label smooth muscle cells in arteriole walls which revealed a blood clot.



In vivo image of mouse ileum microbiome.



Triple labeling of Rat lung. Lung endothelial cell nuclei: acridine orange, 0.05% i.v. Plasma: FITC-Dextran 70kDa, i.v. Distal epithelial membranes: FITC-R.Communis lectin



ViewnVivo enables investigation of living systems in stunning morphological detail.

An advanced tool for longitudinal studies of physiological and pathological processes and cellular, subcellular, and molecular events to answer questions that cannot be visualized by other instruments.

The functional and molecular targeted imaging capability of ViewnVivo® enables the capture of specific cellular events impossible to recreate in-vitro.

Testimonials

“ I had the opportunity over many years to see technologies that come through the laboratory and I would have to say that this probably is the most exciting technology that I have seen in my career come through the laboratory.”

Prof. Mark C Preul, MD

Director of Neurosurgery Research
Barrow Neurological Institute,
Phoenix, Arizona, USA.

“ We have enjoyed using the ViewnVivo system in our labs to study our cells in 3D culture in real time. The system allows us to interrogate the cells and investigate their distribution, proliferation and growth within the collagen matrix at any time in the culture process without removing them from the sterile environment. Using a range of stains we have been able to perform quality control on systems that until now, we were not able to monitor without stopping the experiment and sacrificing samples. The system is well adapted to being integrated into our PC2 hood environment to fit within our standard laboratory process flow.”

Dr Sally McArthur

Director, Manufacturing Futures Research Institute
Swinburne University of Technology,
Melbourne, Australia.

“ The ViewnVivo confocal laser endomicroscope delivers crisp images of ligament and tendon fiber structure simultaneous with fluorescence imaging in a robust and easy to use package. Do you want the ability to collect “big microscope” data where only a millimetre-scale object can fit? Optiscan’s CLE is a powerful tool for tissue research and medical diagnosis.”

Prof. Mark M. Banaszak Holl

Professor and Head
Department of Chemical Engineering
Monash University, Melbourne, Australia.






Image of Glioblastoma stained with fluorescein sodium. The left image shows the margin between tumour and normal brain tissue, while the right image shows the main glioblastoma tumour.



Human Dermal Fibroblasts in collagen gel stained with Alexa 488 Phalloidin (F-actin stain for cell cytoskeleton).



Autofluorescence imaging (left) and confocal endomicroscopy (right) of the anterior cruciate ligament in a study examining failure from fatigue related microdamage.

Technical Specifications

Laser power	10-1000μW
Laser	488nm
Resolution	0.55μm lateral; 5.1μm axial
Field of View	475μm x 475μm
Image Capture Modes	Single frame Continuous capture Z-stack Roll-back (60 frames)
Z-depth	0 – 400μm with dynamic depth change capability
Probes	44mm X 4.0mm Ø 66mm X 4.0mm Ø 300mm X 4.0mm Ø
Filters	8 Standard filters 4 custom filter position
Frame rate	Up to 3.5fps



Footswitch Control



- 1. Home Position and Standby Mode
- 2. Z Depth and Direction Controls
- 3. Image Capture Controls




Precision Micro Positioner



Small probes can be secured in the high precision 3 axis micro positioner mounted to an articulated holding arm for maximum flexibility and accuracy of imaging.



Detection Filters



Band Pass 515-550	Long Pass 515
Band Pass 515-575	Long Pass 540
Band Pass 515-575	Long Pass 570
Band Pass 515-540	

Software Use

Laser Control
Start, pause or adjust the output of the laser

Zoom Control
Provides the ability to change zoom

Scanning Controls
Start or pause scanning



Z Depth Control
Adjustable to navigate through the tissue

Control Panel
Virtual lexan for advanced imaging functions including continuous capture

Brightness Controls
Adjust or activate auto brightness

Filter Selection
Choose from a selection of filters installed

Global Sales Office

Optiscan Imaging Ltd.

T: +61 (0) 9538 3333
E: info@optiscan.com
A: 16 Miles Street, Mulgrave, VIC 3170, Australia

Distributor Information

CHINA, HONG KONG & MACAU



Sinsi Technology Co Ltd.
森西万通科技 (北京) 有限公司

T: +86-10-61666616
E: info@sinsitech.com
A: 27F, Block B, China Railway Construction Plaza, No.19, Beiyuan East Road, Chaoyang District, Beijing, China


SHANGHAI



Biotimes Technology Ltd.
佰泰科技有限公司

T: +86-17701721954 +86-021-54566520
E: liugang@biotimestech.com
A: 24H Aihe Building, No.629 Lingling Road, Xuhui District, Shanghai, China





Optiscan Imaging Limited specialises in Fluorescence In Vivo Endomicroscopy and is the only manufacturer of confocal endomicroscopes with dynamic z-scanning capability.

The flexible hand-held ViewnVivo® probe allows you to push the boundaries of your research with maximum flexibility, as images can be captured at any angle with sub-micron resolution and dynamic optical sectioning.

Global Office 16 Miles Street,
Mulgrave, VIC
3170, Australia

US Office 95 Third Street -
Second Floor, San
Francisco, CA
94103, USA

Email info@optiscan.com

Website



optiscan.com

Optiscan^o
Immediate Informed Decisions